TORAY

HEMOFEEL CH

(English)	Instructions for Use of HEMOFEEL CH Series	P 2~4
(Español)	Instrucciones de uso de HEMOFEEL Serie CH	P 5~7
(Français)	Mode d'emploi de la série HEMOFEEL CH	P 8~10
(Italiano)	Istruzioni per l'uso di HEMOFEEL Serie CH	P 11~13

Instructions for Use of HEMOFEEL CH Series

Read these instructions carefully before using the "HEMOFEEL" CH series haemofilter.

I. DEVICE DESCRIPTION

"HEMOFEEL" CH is a haemofilter, designed for the continuous replacement of the main activity of the kidneys by removing water and toxins at a slow and steady rate (i.e., continuous renal replacement therapy [CRRT]). Haemofilters are devices intended to perform haemofiltration. Haemofiltration is a process whereby extracorporeal blood is slowly and continuously passed through a system to allow for the removal of toxins and/or the replacement of electrolytes. The device functions with dialysate solution for the transfer of substances by diffusion and convection.

II. INTENDED USE

"HEMOFEEL" CH is a single-use haemofilter designed for the continuous replacement of the main activity of the kidneys by removing water and toxins at a slow and steady rate (i.e., continuous renal replacement therapy [CRRT]). The haemofilter should be used only on the direction of a physician.

III. INDICATIONS

"HEMOFEEL" CH haemofilter is used for patients who have acute kidney failure, with conservative treatment consisting of volume, electrolyte and acid-base homeostasis management, and specific drug management.

IV. CONTRAINDICATIONS

Special contraindications for "HEMOFEEL" CH are unknown. Generally, the contraindications for haemofiltration are applicable.

V. STERILIZATION

"HEMOFEEL" CH haemofilters are gamma-ray sterilized and nonpyrogenic (blood side).

VI. PERFORMANCE AND SPECIFICATIONS

Refer to the attached technical data.

VII. ADVERSE EVENTS

- 1) Patients should be carefully monitored during and after haemofiltration. In particular, patients with a history of allergy and hypersensitivity should be carefully monitored. If anaphylactoid reactions should occur, haemofiltration should be discontinued and the attending physician should decide whether or not to return the blood in the extracorporeal circuit to the patient. Symptoms or signs of anaphylactoid reaction includes flushed face, angioedema, urticaria, etc.
- 2) Disequilibrium syndrome may appear especially during haemofiltration treatment at high flow rates. Appropriate treatments, including discontinuation of haemofiltration, should be taken, if the following symptoms or signs of disequilibrium syndrome are observed: headache, nausea, vomiting, consciousness disturbed, decrease of blood pressure or convulsion.
- 3) Appropriate treatments, including discontinuation of haemofiltration, should be taken, if the following symptoms, which may occur occasionally in usual haemofiltration, are observed: Dyspnoea, chest pain, blood pressure decreased, blood pressure increased, thrombocytopenia, headache, nausea, vomiting, abdominal pain, low back pain, leukopenia, palpitation, pyrexia, chills, sweating abnormal, muscle cramp, hypoproteinaemia, hypoalbuminaemia. Patient body fluid and electrolyte levels should be monitored.
- 4) Improper handling or storage during transportation or operation (mechanical shock, excessively high or freezing temperature, etc.) may increase the risk of blood leakage. When blood leakage is found, haemofiltration should be discontinued and appropriate actions should be taken. The attending physician should decide whether or not to return the blood in the extracorporeal circuit to the patient.
- 5) Other complications may occur caused by malfunction, a defect of haemofiltration machine or a procedural error in its use, such as blood loss, blood overheating, hemolysis, excessive ultrafiltration and electrolyte imbalance. Refer to the instructions for the machine to avoid such complications.

VIII. WARNINGS AND PRECAUTIONS

1. General warnings and precautions

- 1) Reuse
"HEMOFEEL" CH is a medical device intended for SINGLE USE ONLY. Its reuse might cause product damages (blood leakage, performance changes, etc.) and also induce patient injury (infection caused by microbiological

contamination, adverse reactions caused by residual reprocessing agents and/or blood components, etc.).

- 2) Rinse
"HEMOFEEL" CH haemofilters are filled with sterile water. Before starting haemofiltration, the blood side of the haemofilters must be flushed with sterile physiological saline, and the dialysate side must be flushed with sterile physiological saline or dialysate to avoid hemolysis, etc.
 - 3) Machine
Since "HEMOFEEL" CH haemofilters have a high ultrafiltration rate, they should be used in conjunction with haemofiltration machines equipped with an ultrafiltration rate controller.
 - 4) Storage
"HEMOFEEL" CH should be stored and transported at temperatures between 5°C (41°F) and 30°C (86°F). Exposure to severe environmental conditions (direct sunlight, excessive humidity, wetness, dryness, freezing temperature) should be avoided.
 - 5) Handling
"HEMOFEEL" CH is made from various plastic components. Excessive vibration or mechanical shock should be avoided during handling and operation. Avoid tapping the haemofilter header with a clamp or other metal object to remove air bubbles or to return the blood.
 - 6) Chemicals
Strictly avoid contact between the device and alcohol (sterilant) or other organic solvents. Such contact may lead to distortions or cracks of the device.
 - 7) Blood line
Some of blood lines have manufacturer's original connectors which do not fit to the dialysis ports of "HEMOFEEL" CH. In that case, use suitable adaptors available in the market.
 - 8) Special care
Since the product has high ultrafiltration performance, start with a low blood flow rate and a low water removal rate in the following patients. [If an inappropriate filtrate flow rate/supplementary liquid flow rate is selected, such symptoms as blood-pressure lowered etc. may occur.]
•Patients with hypotension or in shock
•Patients with low body weight or elderly patients
•Patients in whom a dialyzer with a smaller membrane area was used prior to the use of this product
•Patients who use antihypertensive agents
•Patients who require rapid water removal within a short period
•Patients who experienced blood-pressure lowered etc. due to blood purification therapy
 - 9) Application to elderly patients
Since physiological function generally deteriorates in elderly patients, carefully monitor the patient when using this product.
 - 10) Application to pregnant women or pediatric patients
Since the safety of the product has not been established in pregnant women, use in pregnant women and women of child-bearing potential should be limited to cases where the therapeutic benefits outweigh possible risks.
While the product has been used for pediatric patients, the population is not large enough to establish safety.
- ### 2. Precautions to be taken before haemofiltration
- 1) The haemofilter must not be used if it appears to be damaged or its package is torn.
 - 2) The haemofilter must not be used if there appears to be a water leak such as, drops of water in the package, improper cap fitting over either blood or dialysate port or an accumulation of a large quantity of air in the blood or dialysate compartment.
 - 3) Connect the Luer lock connectors of the blood line to the blood ports by inserting at a correct angle and twisting to ensure a tight fitting to avoid detachment of connection, saline or blood leakage, or breakage of the blood ports.
 - 4) The extracorporeal circuit with the connectors conforming to ISO8638 is recommended for this haemofilter.
 - 5) Connect the dialysate connectors of blood line to the dialysis ports by pushing the line with appropriate force to ensure a tight fitting to avoid detachment of connection, dialysate/filtrate leakage or breakage of the dialysate port.
 - 6) When the connectors of blood line do not fit to the dialysis port of haemofilter, use appropriate commercially available adaptors to connect the blood line to dialysis ports. Ensure tight fittings between connectors of blood line with the adaptor and between the adaptor and the dialysis ports of haemofilter.

7) The haemofilter should be used just after rinsing and priming.

8) Take special care not to allow air bubbles to remain in or enter the blood compartment during the rinsing and priming procedure.

3. Precautions during haemofiltration

- 1) Confirm that there is no blood leakage and no fluid leakage at the blood connection and the dialysate connection sites.
- 2) Transmembrane pressure (TMP) must not exceed 66 kPa (500 mmHg). TMP higher than 66 kPa (500 mmHg) could cause blood leakage or disconnection of the blood lines from the haemofilter connectors.
- 3) The administration method and dose of an anticoagulant such as heparin should be determined for each patient by a physician. Special cares, including monitoring of coagulation time during haemofiltration, should be taken for patients who have a tendency to haemorrhage or who have a coagulation disorder.
- 4) If air bubbles, a blood leak, clotting or hemolysis are observed during haemofiltration, appropriate treatments, including discontinuation of haemofiltration or replacing the haemofilter, should be taken promptly.
- 5) When the administration of drugs is intended before or during haemofiltration, consider various conditions such as route, timing and dose, because haemofiltration may remove them or modify their effects. Special cares are recommended for patients receiving angiotensin-converting enzyme inhibitors.
- 6) When completing haemofiltration, air rinsing of remaining blood should not be employed. Return the patient's blood by rinsing with saline, returning as much of the patient's blood as possible by holding the haemofilter vertically and rotating it around its axis, if necessary.

IX. HAEMOFILTRATION PROCEDURE

Use protective gloves and an aseptic technique to prevent contamination of patient's blood when connecting the blood line to the haemofilter, collecting blood, or returning the blood from the extracorporeal circuit.

The following procedure is an example.

1. Preparation for Haemofiltration

- 1) Introduce physiological saline solution into the end of the arterial side (inlet: Side A in the figure) blood circuit. Avoid air bubbles becoming trapped in the header, and connect the circuit to the arterial side header. Air bubbles in the header must be eliminated by introducing physiological saline solution into the syringe barrel in advance.
- 2) Connect the venous side (outlet: Side V in the figure) header with the venous side blood circuit. (Figure 1)
- 3) Turn the haemofilter upside down so that the arterial side is facing downwards. (Figure 2)
- 4) Use at least 500 mL of physiological saline solution for CH-1.3W, and wash in the order of the arterial side blood circuit, inside of the hollow fibers (blood side), and the venous side blood circuit. (Figure 3) (Use at least 900 mL for CH-1.8W.)
- 5) Turn the haemofilter upside down, remove the stopper of the filter-side haemofilter, and connect the filter side nozzle near the arterial side to the filter pressure monitor. (Figure 4)
- 6) Turn the haemofilter upside down, clamp the venous side with forceps, drain at least 500 mL of physiological saline solution from the blood-side through the haemofilter to remove air bubbles in the filter-side and wash. (Figure 5)
- 7) Turn the haemofilter upside down, and fix it into the holder of the device. (Figure 6)
- 8) Prime the haemofilter and blood lines with more than 500 mL of anticoagulant-added physiological saline, or 5% dextrose solution and confirm that there are no air bubbles remaining in the haemofilter.

2. Haemofiltration

- 1) Connect the blood circuit to the patient's blood removal side catheter. Pay due attention to prevent blood contamination when connecting the blood circuit to the catheter.
- 2) Operate the blood pump slowly; let a small amount of priming solution in the circuit flow and discard the solution while monitoring the patient's condition. When air bubbles in the end of blood return side circuit have been removed, connect the circuit immediately to the patient's blood return side catheter.
- 3) Start haemoperfusion, and add anticoagulant to the blood circulating in the circuit before it enters the haemofilter.

Instructions for Use of HEMOFEEL CH Series

- 4) Maintain a blood flow rate of 20 to 30 mL/min until the entire priming solution in the blood circuit is replaced. When the haemofilter and circuit have become filled with blood, raise the blood flow rate to 50 to 60 mL/min. Perform haemoperfusion alone for at least 5 to 10 minutes after initiation for blank operation without blood filtration.
- 5) After at least 5 to 10 minutes after initiation of haemoperfusion, confirm that haemoperfusion is stable, and then initiate haemofiltration.
- 6) Increase the blood flow rate gradually, and select an appropriate blood flow rate for extracorporeal circulation.
- 7) Infuse anticoagulant continuously to prevent blood from coagulating in the haemofilter while the product is in use. (Figure 7)
- 8) Remain vigilant regarding any difference in pressure between inlet and outlet of the haemofilter – a large difference may be caused by a clogged haemofilter due to blood clot formation etc.

3. Completion of Haemofiltration

- 1) Circulate physiological saline solution or 5% dextrose solution 100 to 200 mL at low flow rate (approx. 50 mL/min) while monitoring the patient's condition, and return the blood in the haemofilter and blood circuit to the patient. If the volume of physiological saline solution used for retransfusion is insufficient, not all the blood may be collected. At the final stage of blood collection, natural pressure is recommended in order to prevent air from becoming trapped in the blood vessel.
- 2) Discard the haemofilter, the lines and other disposables in the manner approved by your institution.

X. WARRANTY AND LIMITATION OF LIABILITY

- 1) "HEMOFEEL" CH is manufactured in accordance with its specifications and in compliance with the GMP regulatory guidelines. Toray will replace a haemofilter with manufacturing defects free of charge. When


necessary, Toray will request return of the haemofilter itself and all related packaging in order to investigate the cause and to take corrective actions, if applicable.

- 2) Toray is not responsible for any damages, irrespective of product damages, patient injury or any other problems caused by misuse, improper handling, operation and storage by the customer.
- 3) Toray is not responsible for any product damages, patient injury or any other problems caused by reuse of the haemofilter.

XI. OTHER INFORMATION

The following information is available on request.

- 1) Details of the test methods
- 2) *In vivo* performance characteristics

						
Manufacturer	Authorized Representative in the European Community	Date of manufacture (Sterilization date)	Use by	Lot No.	Gamma ray sterilized (Sterile fluid path)	Handle with Care
						
Keep away from sunlight	Keep Dry	Temperature limit	For single use only	Caution: Refer to instructions	Non-pyrogenic (Blood side)	

Technical Data

Type	CH-W Series	
	CH-1.3W	CH-1.8W
Housing	Polystyrene	
Fibers	Polymethylmethacrylate (PMMA)	
Inside diameter (μm)	240	
Membrane thickness (μm)	30	
Effective surface area (m ²)	1.3	1.8
Potting material	Polyurethane	
O-ring	Silicone rubber	
Sterilization	Gamma-ray Irradiation	
Blood volume *1 (mL)	94	130
Sieving coefficient *2		
Albumin	0.004	
Inulin	1.0	
β ₂ -microglobulin *3	—	—
UFR <i>in vitro</i> *4 (mL/hr/mmHg)	24	30
Pressure drop of blood compartment*5		
ΔP _B (mmHg)	27.3	19.8
Max. TMP (kPa(mmHg))	66(500)	
Range of blood flow rates (mL/min)	100-300	

*1 Blood volumes are typical data.

*2 Sieving coefficients are typical measured data with bovine blood plasma.
TP:5.5~6.5g/dL, Q_B:200±4mL/min, Q_F:30.0±3.0mL/min/m², Temp:37±1°C


*3 PMMA membrane adsorbs β₂-microglobulin, thus sieving coefficient of β₂-microglobulin is unmeasurable.

*4 UFRs are typical measured data with bovine blood.


Ht: 32±2%, TP:5.5~6.5g/dL, Q_B: 200±4mL/min, Temp: 37±1°C, TMP: 13.3±1.3 kPa (100±10 mmHg)

*5 Pressure drops are typical measured data with bovine blood.

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EC REP Toray International Italy S.r.l.
Toray Medical Co., Ltd.
Toray Industries, Inc.